

Protocol: Bisulfite treatment of DNA, cloning & sequencing to verify specific gene region

DNA Samples:

- Fresh tissue patient samples:
 - ⇒ **Precancerous lesion**
 - ⇒ **Normal colon mucosa**

Gene region:

- *SFRP4*

Amount of DNA for Bisulfite conversion:

Recommended amount: **200-500ng in 20µl** for optimal results.

required concentration: $500\text{ng}/20\mu\text{l} = 25\text{ng}/\mu\text{l}$

required volume: **20µl**

Amount of DNA for PCR:

Template bisulfite converted DNA **~50ng/µl**

1. Bisulfite Treatment of DNA

- 1.1. Ensure DNA samples have thawed completely
- 1.2. Spin down samples and reagents
- 1.3. Add **130µl** of **Lightning Conversion Reagent** to 20µl of DNA sample in a PCR tube. **Vortex & centrifuge** briefly to ensure there are **no** droplets in the cap or sides of tube
- 1.4. Place the tube in a thermal cycler and perform the following steps:
 - I. **98°C for 8 minutes**
 - II. **54°C for 60 minutes**
 - III. Hold at **4°C** for 3-4 minutes and then proceed immediately with purification
- 1.5. Add **600µl** of **M-Binding Buffer** to a *Spin Column* and place the column into a *Collection Tube*
- 1.6. Load the sample from step 2 into the *Spin Column* containing the M-Binding Buffer. Close the cap and mix by inverting several times.
- 1.7. Centrifuge at full speed for **30 seconds**. Discard the flow-through.
- 1.8. Add **100µl** of **M-Wash Buffer** to the column. Centrifuge at full speed for **30 seconds**.
- 1.9. Add **200µl** of **L-Desulphonation Buffer** to the column and let stand at room temperature for **20 minutes**.
- 1.10. After the incubation centrifuge at full speed for **30 seconds**.
- 1.11. Add **200µl** of **M-Wash Buffer** to the column. Centrifuge at full speed for **30 seconds**
- 1.12. Repeat this wash step
- 1.13. Place the column into 1.5ml microcentrifuge tube and add **10µl** of **M-Elution Buffer** directly to the column matrix.
- 1.14. Centrifuge for **30 seconds** at full speed to elute the DNA.
- 1.15. **NanoDrop: Measure Bisulfite Converted DNA concentration to make the dilutions for the PCR reaction - from these calculate the amount of DNA required.**




Figure: Bisulfite conversion of Cytosine to Uracil

2. PCR

- 2.1. Set up PCR machine with conditions (see below)
- 2.2. Calculate amount of each PCR component (see below) for the number of samples
- 2.3. Label the required number of PCR tubes and add the required amount of ddH₂O to each tube
- 2.4. Prepare the master mix of the PCR components
- 2.5. Mix gently by pipetting up and down. **Do not vortex**
- 2.6. Aliquot master mix into sample tubes.
- 2.7. Add **1µl DNA** to the appropriate sample tubes
- 2.8. Briefly spin samples down

Reagent	Volume
PCR master mix	12.5µl
Coral Load	2.5µl
Forward Primer (10µM)	0.5µl
Reverse Primer (10µM)	0.5µl
ddH ₂ O NEG	9µl
ddH ₂ O Sample	8µl

Conditions:

Initial Denaturation: 95°C - 15 min

<i>Denaturation: 94°C - 30 sec</i>] 45 cycles
<i>Annealing: - 30 sec</i>	
<i>Extension: 72°C - 30 sec</i>	

Final Extension: 72°C - 7minutes

4. Gel Electrophoresis

Prepare a **2% gel** in **60mL** volume

- 4.1. Weigh out the required amount of Agarose powder into a conical flask
- 4.2. Fill up to **60ml** with TAE buffer and heat in the microwave until completely melted
- 4.3. Whilst melting, set-up the chamber with the 8 well-comb
- 4.4. Allow to cool by running cold water over the flask
- 4.5. Add Gel Red to the conical flask containing the melted gel
- 4.6. Pour into the chamber and remove any bubbles from the gel
- 4.7. Leave to set
- 4.8. Set up the chamber for the electrodes and fill up with TAE Buffer
- 4.9. Once set, place the gel into the chamber containing the TAE buffer and pour more buffer over the top to ensure the wells are completely submerged
- 4.10. Load the gel with PCR products
- 4.11. Attach the electrodes and allow the gel to run
- 4.12. Visualise gel (picture)

5. TOPO TA Cloning reaction

Reagent	Volume
PCR product	2µl
Salt Solution	1µl
H2O	2µl
TOPO vector	1µl
Final Volume	6µl

- 5.1. Remove Salt solution, H2O and vector from the -20C to thaw on ice
- 5.2. Add the reagents in the order above. Mix gently and incubate for *5 minute* at **RT**
- 5.3. Leave the reaction on ice




Figure: TOPO TA cloning reaction

6. Preparation for Transformation

Components

S.O.C medium
TOPO Cloning reaction from step 4.1
LB plates containing 50 µg/mL ampicillin
42°C water bath
37°C shaking and non-shaking incubator
General microbiology supplies (e.g. plates, spreaders)

- 6.1. Equilibrate water bath or heat block to **42°C**
- 6.2. Warm the vial of S.O.C medium to **RT**
- 6.3. Warm selective plates at **37°C** for *30 minutes*
- 6.4. Thaw on ice 1 vial of E.coli TOP10 chemically competent cells for each transformation

7. Transformation

- 7.1. Add **2µl** of the TOPO Cloning reaction from **step 4.1** into a eppendorf tube of the chemically competent E.coli and mix gently. (Do not mix by pipetting up and down, instead gentle flicking of the tube)
- 7.2. Incubate on ice for *30 minutes*
- 7.3. Heat shock the cells for *30 seconds* at **42°C** without shaking
- 7.4. Immediately transfer the tubes to ice
- 7.5. Add **250µl** of room temperature S.O.C medium
- 7.6. Cap the tube tightly and shake the eppendorf tube horizontally (200rpm) at **37°C** for *1 hour*
- 7.7. Spread **10 - 50µl** from each transformation on a pre-warmed selective plate and incubate overnight at **37°C**

8. Analysis of Transformants

- 8.1. Prepare a 14ml tube (with a vent cap) for each colony with 4ml of LB medium containing 100µg/mL ampicillin
- 8.2. Pick single colonies with a pipette tip and put into tubes
- 8.3. Incubate overnight at **37°C, 250rpm**

9. Isolation of the plasmid DNA using Plasmid Mini-Prep kit

- 9.1. Centrifuge (at **RT, 5000rpm, for 10 minutes**) the 14ml tubes containing **4ml** of the overnight LB-culture
- 9.2. Remove all of the medium: Pour the majority into a flask for autoclaving and remove the rest with a pipette.
- 9.3. Pre-heat an aliquot of TE buffer to **65-70°C**

- 9.4. Add **250µl** *Resuspension Buffer* (R3) with RNase A to the cell pellet and re-suspend the pellet using a pipette until it is homogenous. Transfer to eppendorf tubes.
- 9.5. Add **250µl** *Lysis Buffer* (L7). Mix gently by inverting the capped tube until the mixture is homogenous. DO NOT vortex. Incubate the tube at **RT** for *5 minutes*.
- 9.6. Add **350µl** *Precipitation Buffer* (N4). Mix immediately by inverting the tube, or for large pellets, vigorously shaking the tube, until the mixture is homogenous. DO NOT vortex. Centrifuge the lysate at $>12,000 \times g$ for *10 minutes*.
- 9.7. Load the supernatant onto a spin column in a 2mL wash tube. Centrifuge the column at $12,000 \times g$ for 1 minute. Discard the flow-through and place the column back into the wash tube.
- 9.8. Add **700µl** *Wash Buffer* (W9) with ethanol to the column. Centrifuge the column at $12,000 \times g$ for *1 minute*. Discard the flow-through and place the column into the wash tube. Centrifuge the column at $12,000 \times g$ for *1 minute*. Discard the wash tube with the flow-through.
- 9.9. Place the spin column in a clean 1.5ml eppendorf. Add **75µl** of preheated *TE buffer* (TE) to the centre of the column. Incubate the column for 1 minute at RT.
- 9.10. Centrifuge the column at $12,000 \times g$ for *2 minutes*. The Eppendorf contains the purified plasmid DNA. Discard the column. Store the plasmid DNA at 4°C (short-term) -20°C (long term)
- 9.11. Measure concentrations on the NanoDrop

10. Analyse the transformants by restriction enzyme digest

- 10.1. Make a master mix for the following for the number of samples required:
 - 1.5µl of NEBuffer 2.1
 - 0.2µl of EcoR1
 - 12.3µl of H₂O
- 10.2. Add 14µl of the master mix for the number of samples
- 10.3. Add 200ng of DNA to the samples with master mix
- 10.4. Incubate for at least 1 hour at 37°C
- 10.5. Prepare gel
- 10.6. Run transformants on gel and visualise

11. Sequencing

11.1. Prepare samples to be sent to Microsynth for sequencing

